

Instrument Set 3

Micro Flap Management Set



Secure and lasting wound closure is an essential factor for successful interventions in augmentation surgery.

The success of the therapy greatly depends on the design and mobilization of the flap used for covering the wound, and no less important is precise microsurgical wound edge adaptation to ensure the best possible blood supply.

Coming with a whole range of high-quality microinstruments, the micro flap management set provides everything needed for precise and successful soft tissue management.



Set 3: Micro Flap Management Set
consisting of:

Item No.	Qty	Designation
10-160-01-07	1	Scalpel handle for microblades, 13.5 cm
10-120-04-07	1	Scalpel handle no. 3, 16 cm
23-402-18-07	1	Periosteal elevator, micro, Iglhaut, 18 cm
43-300-01-07	1	Periodontal probe, Iglhaut
20-100-18-07	1	Micro needleholder, Iglhaut, straight, with ratchet, 18 cm
12-562-18-07	1	Micro forceps, Iglhaut, 1 x 2 teeth, straight, 0.7 mm, 18 cm
12-563-18-07	1	Micro forceps, Iglhaut, 1 x 2 teeth, curved, 0.7 mm, 18 cm
12-560-18-07	1	Micro forceps, Iglhaut, straight, 6.0 x 0.7 mm, 18 cm
11-715-18-07	1	Micro scissors, Iglhaut, curved

MiniSet sterile container MicroStop®		
55-861-50-04	1	Container 310 x 189 x 70 mm
Accessories container		
55-864-01-04	1	Coding label for cover, with text: Micro flap management set
55-864-05-04	1	Coding label for front, with text: Micro flap management set
55-864-15-04	1	Logistics framelets for cover, yellow

Mesh tray		
55-015-10-01	1	Mesh tray 277 x 171 x 30 mm
Accessories mesh tray		
55-004-03-04	2	Bottom strip for instrument fixation, transverse, 162 mm
55-005-01-04	2	Plug-in element with support (16/pack)
55-006-02-04	1	Spring element, transverse, 162 mm
55-008-01-04	1	Stop frame, transverse, 162 mm








Fig. 10 - 17

10-160-01-07
13.5 cm / 5 3/8"
Scalpel handle for microblades

10-120-04-07
16.5 cm / 6 2/8"
Scalpel handle no. 3



Iglhaut
23-402-18-07
18 cm / 7 1/8"
Micro periosteal elevator



Iglhaut
43-300-01-07
17.5 cm / 6 7/8"
Periodontal probe,
solid handle, 6 mm



1/2



1/1



diamond-coated



safety-wave



4-0 - 7-0

Igihaut
20-100-18-07
18 cm / 7 1/8"
Micro needleholder



1/2



1/1

12-560-18-07
0.70 mm



1/1

12-562-18-07
0.70 mm



1/1

12-563-18-07
0.70 mm



1:2



1:2

Igihaut
12-560-18-07 - 12-563-18-07
18 cm / 7 1/8"
Micro forceps




1/2



1/1



10 mm



safety-wave

Igihaut
11-715-18-07
18 cm / 7 1/8"
Micro scissors