



SINGLE USE

SUPER IMAGE QUALITY

VARIETY OF MODELS

Video processor incorporates both freeze & video capture to USB



Imaging Processor
PL-2100(Mini)



Imaging Processor
PL-2000

SPECIFICATIONS

	PL-1000	PL-2000 (Portable)	PL-2100 (Mini)
Power	Voltage 100-240 VAC Frequency 50/60 Hz Rated input 45 VA	Voltage 100-240 VAC Frequency 50/60 Hz Rated input 24 VA	Voltage 100-240 VAC Frequency 50/60 Hz Rated input 24 VA
Illumination	Lamp 2W LED Lamp life 30,000 hrs Color temperature 3,000K-7,000K	Lamp 2W LED Lamp life 30,000 hrs	Lamp 2W LED Lamp life 30,000 hrs
Function	White balance, Luminance control(+/-), Sharpness(0.1,2), Zoom(0.1,2), Freeze, Snap, Record, USB storage	White balance, Luminance control(+/-), Sharpness(0-10), Zoom(0.1,2), Freeze, Snap, Record, USB storage	White balance, Luminance control(+/-), Sharpness(0-10), Zoom(0.1,2), Freeze, Snap, Record, USB storage
Output	DVI*1, S-Video*1, CVBS*1	HDMI*2	HDMI*2
Storage	USB 2.0 (FAT32, Max.128G)	USB 2.0 (FAT32, Max.128G)	USB 2.0 (FAT32, Max.128G)
Dimensions	350mm(W) * 130mm(H) * 425mm(D)	390mm(W) * 270mm(H) * 200mm(D)	270mm(W) * 260mm(H) * 72mm(D)
Weight	8 kg	7.2 kg	3 kg
Screen	/	12.1 inch	10.1 inch

Imaging Processor
PL-1000



Distributed by



Shanghai SeeGen Photoelectric Technology Co.,Ltd.

Add: 3 Floor, Building No.1 4299 Jindu Road,
Minhang Dist., Shanghai China

www.seegenmed.com

info@seegen.com.cn


+86 021 5429 5355

「Briview」


Disposable Choledochoscope
- for ERCP use

Shanghai SeeGen Photoelectric Technology Co., Ltd.

Choledochoscope - for ERCP use




SUPERIOR
High resolution image



Higher resolution image.
Twin light guides deliver clear, bright image.

4 way deflection
on treatment scope.

Zoom
Freeze



Working with duodenoscope

	Diagnose	Treatment
Mode	CHV-US100J-U	CHV-US110J-U
Field of view	110°	110°
Depth of view	5-50mm	5-50mm
Deflection	up/down 120°/120°	up/down 120°/120° left/right 120°/120°
Insertion dia.	3.3mm	3.4mm
Distal end dia.	3.3mm	3.4mm
Working length	1870mm	1870mm
Working channel	1.2mm/0.6mm	1.2mm/0.6mm

Water jet



Suction

Working length: **1870mm**