

TWISTER® PLUS

Rotatable Retrieval Device



Horizons International Corp.

P.O. Box 213-3006 Zona Franca Metro, Barreal de Heredia, Costa Rica
Phone (506) 2261-9114 • Fax (506) 2261-9115
www.horizonscorp.com • info@horizonscorp.com


TWISTER® PLUS

Rotatable Retrieval Device

TWISTER® PLUS Rotatable Device is a fully rotatable, three dimensional net that is designed to facilitate polyp, food and foreign body capture and removal.



FEATURES AND BENEFITS



360° ROTATION (1:1 RATIO IN ALL DIRECTIONS)

Designed to facilitate capture and removal of polyps, food boluses and foreign bodies

DEEP, THREE- DIMENSIONAL NET

Engineered to capture and remove multiple polyps, large food boluses and foreign bodies in a single pass

DURABILITY: INTERWOVEN NET DESIGN

Net is attached to the outer wire using a interwoven design to help maintain the integrity of the net

SHORT THROW HANDLE

Ergonomically design for nurse and technician comfort

FOUR NET SIZES

15 mm, 22 mm, 26mm and 32 mm net sizes are available to meet specific procedural requirements

PRODUCT SPECIFICATIONS

ORDERING INFORMATION						
Product Code	Description	Net Size (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	Packing
HIC-GDL-27-5	TWISTER® PLUS Rotatable Retrieval Device	22	2.5	230	2.8	Box 5
HIC-GDL-27-10	TWISTER® PLUS Rotatable Retrieval Device	22	2.5	230	2.8	Box 10
HIC-GDL-26-5	TWISTER® PLUS Rotatable Retrieval Device	26	2.5	230	2.8	Box 5
HIC-GDL-26-10	TWISTER® PLUS Rotatable Retrieval Device	26	2.5	230	2.8	Box 10
HIC-GDL-29-5	TWISTER® PLUS Rotatable Retrieval Device	32	2.5	230	2.8	Box 5
HIC-GDL-29-10	TWISTER® PLUS Rotatable Retrieval Device	32	2.5	230	2.8	Box 10
HIC-GDL-30-5	TWISTER® PLUS Rotatable Retrieval Device	15	1.7	120	2.0	Box 5
HIC-GDL-30-10	TWISTER® PLUS Rotatable Retrieval Device	15	1.7	120	2.0	Box 10
HIC-GDL-31-5	TWISTER® PLUS Rotatable Retrieval Device	15	1.7	160	2.0	Box 5
HIC-GDL-31-10	TWISTER® PLUS Rotatable Retrieval Device	15	1.7	160	2.0	Box 10

TWISTER® PLUS is a registered Trademark of Horizons International Corp.
Protected by one or more of the following US Patents: 9,101,342; 8,858,567.

Caution: Federal (USA) law restricts this device to be sold by, or on the order of a physician.