



Pelican
feminine healthcare
an eakin company 

Pages	Contents
10	How To Order
11-13	NHS Order Codes
14-20	PELIspec® Vaginal Speculum Range
21-23	Proctoscope & Lighting System
24-25	Gynaecological & Instrument Lubricant
26-30	Scissors & Forceps
31-35	Intra-Uterine Device (IUD) Removal & Fitting
36-39	Cervical Dilator & Uterine Sound
40-42	Obstetrics Range
43-47	Model-Med: Anatomical Training Models

Welcome to Pelican Feminine Healthcare

Pelican Feminine Healthcare Ltd is a UK manufacturer of **single-use** medical instruments for **Gynaecology** and **Obstetrics**, including the market leading PELIspec® vaginal speculum.

State-of-the-art manufacturing facilities and active research and development highlight Pelican Feminine Healthcare's commitment to the continuous development of **innovative healthcare solutions** in Women's health.

This, combined with a positive input from specialist clinicians and other professional customers, ensures that the company remains at the forefront of providing **clinically and cost effective products.**

Pelican
feminine healthcare



Innovative Healthcare Solutions

Delivering Quality, Value & Trust



Pelican Feminine Healthcare are at the forefront of medical device innovation in Gynaecology and Obstetrics. Through attention to detail, materials sourcing and exceptional manufacturing capabilities, devices such as the PELIspec® are the market-leading products for health professionals around the world. Pelican FH continues its commitment to innovative products, with the FETOfit® and PELIspec® Pro-Wall, engaging with health professionals to identify unmet clinical needs and develop a new product's pipeline that will add to the high quality, trusted devices that Pelican are recognised for.

Pelican FH partner with NHS Trusts, Academic Partnerships, Universities, Start-ups and Individuals to Research products that make a real difference to the professionals who use them and the people who receive the care required. Through taking proactive measures to ensure that products are BPA and Latex-free, Pelican FH set themselves apart and ensure that quality is at the heart of what we do.

Pelican FH's innovative products have revolutionised the single-use market and are the products of choice for health services and professionals worldwide.

Through innovation Pelican FH have brought products that reflect the quality and value required by clinicians and ensures that they can trust the device in their hand. We act with integrity in everything that we do and provide a service that is second to none.



15,000,000
Market-leading
PELIspec[®]
vaginal speculum
manufactured
since 2011

Pelican Feminine Healthcare
Production Statistics 2011 to 2014

Quality Assurance

Committed to providing quality products, manufactured to the highest standards.

Delivering market leading innovations to the highest quality standards is reinforced by **stringent monitoring** of manufacturing processes with over 15,000 quality control checks performed annually and across Pelican Feminine Healthcare's product range.

Pelican FH has invested heavily in the design and manufacture of products including the market-leading PELIspec® - **strength tested every four hours** to ensure it meets stringent product specifications. These checks are conducted at a state-of-the-art manufacturing facility in Cardiff, which since 2011 has served both UK and international markets with production of over fifteen million vaginal speculum.

The company is focused on providing quality products, manufactured to the highest standards, using **superior medical grade materials** and continuous improvement is the hallmark of Pelican FH's activities.



Pelican Feminine Healthcare are full accredited for ISO 13485:2003 and all products are CE marked in accordance with European Directive 93/42/EEC. All sterile processes are assessed for sterility by an MHRA accredited regulatory body.

Each device has a unique product identification number which can be recorded on patient notes for traceability.

All products are Latex-free, sterile, and comply with National Decontamination Guidelines.

15,000+

Quality control
checks performed
annually across
the product range

Pelican Feminine Healthcare
Production Statistics 2011 to 2014





How To Order

Browse Online:

Visit www.pelicanfh.co.uk to explore our latest range of products.

Via NHS Supply Chain:

Our most popular products are available through NHS Supply Chain. Please see Page 11-13 for details of applicable NHS Order Codes.

Via UK Distributor Partners:

For more information on UK regional distributors of Pelican Feminine Healthcare products, please contact us: contactus@pelicanfh.co.uk, Freephone: 0800 013 1220.

**Order Direct**

Freephone: 0800 013 1220

Fax: 02920 767 801

E-mail: telesales@pelicanfh.co.uk

**International Orders**

Please call us: +44 (0) 2920 767 800

or email: export@pelicanfh.co.uk

Delivery

Delivery is to mainland England and Wales only. Please note that we do not currently deliver to Scotland or Southern Ireland please contact us on Freephone 0800 013 1220 for our distributor partners in these areas.

Orders placed **before 3pm** are **guaranteed next day delivery** within the UK and two working days for Northern Ireland.

NHS order codes

Our specialist range of feminine healthcare products are assigned NHS Order codes so can be easily ordered directly from NHS Supply-Chain.



Code	Description	Pack Size	NHS Supply Chain Order Code
400091	PELIspec® Vaginal Speculum - Virgin with Lock	25 Units	FGR023
400100	PELIspec® Vaginal Speculum - Small without Lock	25 Units	FGR447
400101	PELIspec® Vaginal Speculum - Medium without Lock	25 Units	FGR022
400102	PELIspec® Vaginal Speculum - Medium Long without Lock	25 Units	FGR448
400103	PELIspec® Vaginal Speculum - Mixed with Lock	25 Units	FGR089
400104	PELIspec® Vaginal Speculum - Mixed without Lock	25 Units	FGR449
400105	PELIspec® Vaginal Speculum - Small with Lock	25 Units	FGR019
400106	PELIspec® Vaginal Speculum - Medium with Lock	25 Units	FGR013
400107	PELIspec® Vaginal Speculum - Medium Long with Lock	25 Units	FGR020
400108	PELIspec® Vaginal Speculum - Broad with Lock	25 Units	FGR008
400111	PELIspec® with Smoke Extractor - Medium	20 Units	FGR002
400112	PELIspec® with Smoke Extractor - Medium Long	20 Units	FGR005
400113	PELIspec® with Smoke Extractor - Broad	20 Units	FGR088
400131	PELIspec® Sims Speculum - Medium	25 Units	FGR021
400135	PELIspec® Sims Speculum - Large	10 Units	FGR435
400200	PELIspec® Vaginal Speculum with Light Source (Starter Pack)	1 Unit	FGR431
400207	PELIspec® Vaginal Speculum (Medium) with Light Source	25 Units	FGR432
400209	PELIspec® Vaginal Speculum (Medium Long) with Light Source	25 Units	FGR433
400210	Leads for Light Source Vaginal Speculum	5 Units	FGR437
400220	Battery Pack for Light Source Vaginal Speculum	1 Unit	FGR438
400233	PELIspec® Pro-Wall Speculum Medium	10 Units	FGR472
400235	PELIspec® Pro-Wall Speculum Medium Long	10 Units	FGR473

Code	Description	Pack Size	NHS Supply Chain Order Code
400300	Pelican Spongeholder	50 Units	FGR221
400310	Pelican Uterine Sound	100 Units	FGR030
400317	Pelican Coil Remover	50 Units	FGR035
400320	Pelican Long Scissors	50 Units	FGR032
400330	Pelican Tenaculum	50 Units	FGR029
400335	Pelican Vulsellum	50 Units	FGR006
400342	Pelican Dilator 3/4mm	100 Units	FGR430
400344	Pelican Dilator 5/6mm	100 Units	FGR031
400347	Pelican Dilator 7/8mm	100 Units	FGR434
400351	PELlpack Removal & Fitting Kit - Medium	4 Units	FGR036
400352	PELlpack Removal & Fitting Kit - Medium Long	4 Units	FGR037
400400	PELljelly Lubricant	100 Units	FGR090
400560	Pelican Thread Retriever	5 Units	FGR220
400565	Pelican Thread Retriever	25 Units	FGR436
400625	FETOfit®	25 Units	FGR450
400706	Pelican IUD Pack (Medium)	4 Units	FGR451
400707	Pelican IUD Pack (Medium Long)	4 Units	FGR452
400708	Pelican IUD Fitting Pack (Medium)	4 Units	FGR429
401000	Uterine Sound-Dilator	50 Units	FGR471
KH2	PELlhook Amniotic Membrane Perforator - Curved	100 Units	FFY059
SH2	PELlhook Amniotic Membrane Perforator - Straight	100 Units	FFY060

PELIspec® Vaginal Speculum Range



Fail
Safe



The marketing-leading PELIspec® vaginal speculum range continues to lead the way in providing a high quality single-use speculum to the healthcare market, and has been the speculum of choice for discerning medical professionals for **over 20 years**.

Engineered and moulded in the UK using specially chosen medical grade BPA-free (Bisphenol-A) plastic, the PELIspec® vaginal speculum range provides the ultimate combination of strength and flexibility with patient safety.

“Fail-SAFE” points are uniquely incorporated into the PELIspec® design, ensuring that in the event of ‘Vaginismus’ or other unanticipated stress the device will flex and break safely outside the vagina.

The PELIspec® range offers a **wide variety** of speculum solutions tailored to each procedure's requirements from built-in smoke extraction and light source to the ‘Pro-Wall’ with additional wall support. The more conventional core range accommodates for every patient, available in a wide range of sizes from virgin and medium to medium-long and broad, ‘with’ and ‘without lock’.

Pelican Feminine Healthcare is committed to promoting awareness of cervical cancer and continues to donate **5p from every box of PELIspec® vaginal speculum sold**, to Jo's Cervical Cancer Trust - the only UK charity dedicated to women and their families affected by cervical cancer and cervical abnormalities. We have now **donated over £40,000*** to support this vital cause.

*Figure accurate as of March 2015


Pelican are contributing 5p for every box of PELIspec sold




Registered charity in England and Wales, number: 1133542
Registered charity in Scotland, number: SC041236

www.pelicanfh.co.uk

“Strong, but flexible
design engineered
with safety
in mind”





PELIspec® Vaginal Speculum

The market leading PELIspec® is a strong, sterile, disposable vaginal speculum, available in 5 sizes. Selected sizes are also available without the locking mechanism.

PELIspec® can be used for gynaecological procedures, including cervical screening, IUD fittings and vaginal examinations.

Features

- ✓ Proven strength with built-in 'Fail-Safe' points
- ✓ Textured handles for enhanced grip
- ✓ Crystal clear for maximum visibility
- ✓ BPA and Latex-free
- ✓ Unique product identification number for traceability




Fail
Safe



Manufactured
in the UK



Code	Description	Units
400091	PELIspec® - Virgin with Lock	25
400105	PELIspec® - Small with Lock	25
400106	PELIspec® - Medium with Lock	25
400107	PELIspec® - Medium Long with Lock	25
400108	PELIspec® - Broad with Lock	25
400103	PELIspec® - Mixed Box with Lock	25
400100	PELIspec® - Small without Lock	25
400101	PELIspec® - Medium without Lock	25
400102	PELIspec® - Medium Long without Lock	25
400104	PELIspec® - Mixed Box without Lock	25



PELIspec® Pro-Wall

NEW

The PELIspec® Pro-Wall is the ‘Professional’ solution specially designed for use in cases of anterior vaginal wall laxity.

Based on the market leading PELIspec® vaginal speculum, the device features a Latex-free sheath configured to efficiently secure an unobstructed view of the cervical anatomy to ensure maximum visibility during the procedure.

Features

- ✓ Ready applied Latex-free sheath
- ✓ Proven strength with built-in ‘Fail-Safe’ points
- ✓ BPA and Latex-free
- ✓ Textured handles for enhanced grip
- ✓ Unique product identification number for traceability



Fail
Safe



Manufactured
in the UK

Code	Description	Units
400233	PELIspec® Pro-Wall Speculum - Medium	10
400235	PELIspec® Pro-Wall Speculum - Medium Long	10



PELIspec® with Light Source

Based on the market leading PELIspec® vaginal speculum, the unique PELIspec® with Light Source features a built-in light source attachment configured to provide effective directional light for maximum visibility during gynaecological procedures.

Available in a 'Starter pack' - containing 2 x PELIspec® (1 x Medium, 1 x Medium Long), 2 x Reusable Light source leads, 2 x plastic lead covers, 1 x Re-usable Battery Pack.

Features

- ✓ Directional light source for maximum visibility
- ✓ Compact re-usable battery pack with belt clip
- ✓ Sealed light source insertion point to ensure protection during use
- ✓ Proven strength, with 'Fail-SAFE' Points uniquely built into each PELIspec®
- ✓ Available in two sizes – Medium and Medium Long
- ✓ Reusable Leads supplied with plastic covers to reduce cross-contamination



Fail
Safe



Manufactured
in the UK

Code	Description	Units
400207	PELIspec® with Light Source - Medium*	25
400209	PELIspec® with Light Source - Medium Long*	25
400200	PELIspec® with Light Source Starter Pack - Mixed	1
400210	Light Source Leads**	5
400220	Light Source Battery Pack	1

*Lead covers included/ **Lead covers not included



PELIspec® with Smoke Extractor

The PELIspec® with Smoke Extractor is specially designed for use during loop diathermy procedures.

Based on the market leading PELIspec® vaginal speculum, the device features a built-in extraction tube configured to efficiently remove smoke produced during the procedure.

The PELIspec® with Smoke Extractor is deliberately manufactured in white to best reflect light and ensure maximum visibility.

Features

- ✓ Effective smoke extraction
- ✓ Proven strength with built-in 'Fail-Safe' points
- ✓ BPA and Latex-free
- ✓ Textured handles for improved grip



Manufactured in the UK



Code	Description	Units
400111	PELIspec® with Smoke Extractor - Medium	20
400112	PELIspec® with Smoke Extractor - Medium Long	20
400113	PELIspec® with Smoke Extractor - Broad	20



PELIspec® Sims Speculum

The PELIspec® Sims Speculum is commonly used for inspection of the vaginal walls and can be used as an alternative to the traditional Cusco design.

The design and shape of the PELIspec® Sims Speculum is characterised by its two right-angled bends and differing sized ends.

Features

- ✓ Crystal clear ensuring greater visibility of the vaginal wall
- ✓ Smooth rounded ends for patient comfort
- ✓ Mirrors design of metal equivalent
- ✓ Easy central grip handle
- ✓ Two differing end sizes for convenience



Manufactured
in the UK

Code	Description	Units
400131	PELIspec® Sims Speculum - Medium	25
400132	PELIspec® Sims Speculum - Medium	3
400135	PELIspec® Sims Speculum - Large	10

Proctoscope & Lighting System





Pelican Proctoscope

The Pelican Proctoscope is specially designed for use during proctoscopy and features a crystal clear barrel with a smooth moulded introducer for gentle insertion into the rectum. It is recommended that PELIjelly Lubricant is used in conjunction with the device to minimise discomfort upon insertion.

The Pelican Proctoscope incorporates a housing to allow connection to an external light source which is sealed to minimise the risk of cross-infection.

Features

- ✓ Sealed housing to allow connection to an external light source
- ✓ Crystal clear
- ✓ Latex and BPA-free
- ✓ Unique product identification number for traceability



Manufactured
in the UK



Code	Description	Units
400370	Pelican Proctoscope - Paediatric	25
400380	Pelican Proctoscope - Medium	25
400390	Pelican Proctoscope - Large	25

Pelican Lighting System



Complete fibre optic lighting system for use in conjunction with the Pelican Proctoscope.

Range of adaptors is also available to allow connection to an existing light source.

Features

- ✓ 1.8 metre fibre optic illumination cable
- ✓ 150 watt halogen light source
- ✓ Output intensity adjustment
- ✓ Detachable adaptor on cable to allow sterilisation
- ✓ 12 months warranty
- ✓ Conforms to European standards



Use With Our Proctoscope

Code	Description	Units
400500	Pelican Fibre Optic Lighting System	1
400510	Naughton Morgan/ Storz Olympus Adaptor	1
400520	Heine Adaptor	1
400530	Woolf Adaptor	1
400540	Welch Allyn Adaptor	1
400550	Sapimed Adaptor & Cable	1

Gynaecological & Instrument Lubricant

PELLjelly is a sterile lubricant, ideally suited for gynaecological examinations.

- ✓ **Colourless, odourless and water-based**
- ✓ **Easy to apply, safe to use.**

An article published by Margie Givens, Senior Biomedical Scientist at Forth Park Hospital, examined the use of lubricating gels during cervical screening. A summary of this article concluded that:

“PELLjelly may be better used (than KY Jelly) when taking cervical samples”*

“PELLjelly still had sufficient (cervical) cells for diagnosis. The KY sample did not...”*

The full article is available to download online via the Pelican Feminine Healthcare website

*Givens, M (2008) “Cytology Uncovered! Lubricating Gels - a Personal Experience”. *Obs Gynae & Midwifery, News, Autumn 2008*, pp.16-18.

“Easy to apply,
safe to use”.





PELIjelly Lubricant

PELIjelly is a sterile, water based lubricant designed for use in gynaecological, digital and instrument examinations with no effect on DNA.*

The sachet format allows minimal waste and also removes any concerns regarding cross-infection.

Features

- ✓ Convenient Sterile 5g sachet
- ✓ No effect on DNA*
- ✓ Silk like consistency to minimise discomfort from procedure
- ✓ Colourless and fragrance free
- ✓ Easy to apply and will not stain
- ✓ 100 pack dispensing carton



Manufactured in the UK



Code	Description	Units
400400	PELIjelly Lubricant	100

*Givens, M (2008) "Cytology Uncovered! Lubricating Gels - a Personal Experience". Obs Gynae & Midwifery, News, Autumn 2008, pp.16-18.

Scissors & Forceps





Pelican Spongeholder

The Pelican Spongeholder is designed to ensure visibility is not impeded when used in conjunction with a PELIspec® vaginal speculum.

The robust design gives the user full control of the device and the grooved tips ensure a secure grip of swabs.

Features

- ✓ Grooved tip for maximum grip
- ✓ Enhanced secure locking mechanism
- ✓ Specifically designed for use with swabs
- ✓ Cost effective alternative to metal with parallel performance



Code	Description	Units
400300	Pelican Spongeholder	50
400302	Pelican Spongeholder	5



Pelican Vulsellum

The Pelican Vulsellum is characterised by a 3:4 toothed grip specially designed to offer a gentle but secure grip on the cervix.

The Pelican Vulsellum has an angled head to maximise access to the cervix during the procedure and improve visibility.

Provides an alternative choice for practitioners who prefer a more evenly spread grip over a singled toothed (Tenaculum) grip.

Features

- ✓ Enhanced secure locking mechanism
- ✓ Specifically designed tooth combination to minimise discomfort and reduce trauma
- ✓ Angled tips allow greater visibility
- ✓ Cost effective alternative to metal with parallel performance



Manufactured
in the UK

Code	Description	Units
400335	Pelican Vulsellum	50
400336	Pelican Vulsellum	5



Pelican Tenaculum

The Pelican Tenaculum features a precise single toothed grip with secure locking mechanism, ideal for holding the anterior lip of the cervix during procedures.

The Tenaculum provides an alternative choice to practitioners who prefer a single toothed grip over a more evenly spread grip.

Features

- ✓ Precise single toothed grip
- ✓ Secure locking mechanism
- ✓ Cost effective alternative to metal with parallel performance
- ✓ Specifically designed to allow greater visibility of the cervix



Manufactured
in the UK

Code	Description	Units
400330	Pelican Tenaculum	50
400331	Pelican Tenaculum	5



Pelican Long Scissors

The Pelican Long Scissors are ideal for cutting IUD threads following the introduction of an intra-uterine device (IUD). The slender, flexible handles allow full manoeuvrability of the device, and the blade design ensures that the threads are captured securely for a clean cut.

Features

- ✓ Specially designed for cutting coil threads
- ✓ Sharp blade for a clean cut
- ✓ Firm finger grip for greater control
- ✓ Long handle for ease of use



Manufactured
in the UK

Code	Description	Units
400320	Pelican Long Scissors	50
400321	Pelican Long Scissors	5

Intra-Uterine Device (IUD) Removal & Fitting





Pelican PELIpack - IUD Removal & Fitting Pack

The Pelican PELIpack contains all the instrumentation necessary to remove and fit an intra-uterine device. The instruments are presented in an individual sterile pouch offering a convenient time-saving solution.

All instruments from the Pelican PELIpack are also available separately.

Contents

- PELIspec® Medium (Pack code 400351) or PELIspec® Medium Long (Pack code: 400352)
- Pelican Spongeholder
- Pelican Vulsellum
- Pelican 3/4mm Dilator
- Pelican Uterine Sound
- Pelican Long Scissors
- Pelican Coil Remover
- Gallipot, Gauze, Drape & Tray



Manufactured
in the UK



Code	Description	Units
400351	PELIpack IUD Removal & Fitting Kit - Medium	4
400352	PELIpack IUD Removal & Fitting Kit - Medium Long	4



Pelican IUD Pack - Removal & Fitting Pack

The Pelican IUD Pack is a single-use procedure pack containing a mixture of single-use metal and plastic instruments to remove and fit an intra-uterine device.

Each sterile peel pouch incorporates a peel-off label which can be attached to patient notes.

Contents

- PELIspec® Medium or PELIspec® Medium Long
- Metal Rampley Spongeholder Forceps*
- Metal Teale Vulsellum Forceps
- Metal Hegar Dilator*
- Metal Mayo Scissors
- Pelican Uterine Sound
- Pelican Coil Remover*
- Gallipot*, Cotton Wool Balls*, Drape & Kidney Dish



Manufactured
in the UK

Code	Description	Units
400706	Pelican IUD Pack (Medium)	4
400707	Pelican IUD Pack (Medium Long)	4
400708	Pelican IUD Fitting Only Pack (Medium)	4

*Pack codes: 400706 and 400707 only

Pelican Coil Remover

The Pelican Coil Remover is specifically designed to remove intra-uterine devices and is characterised by an angled tip with a 2:1 tooth combination.

The robust design gives the user full control of the instrument and the “V” designed tip ensures a secure grasp of the coil thread to facilitate retrieval of the intra-uterine device. The tip size has been reduced and angled for optimal visibility during the procedure.

Features


- ✓ Secure locking mechanism
- ✓ Angled tip to optimise visibility
- ✓ Specially designed tooth combination facilitates easy IUD retrieval



Manufactured
in the UK



Code	Description	Units
400317	Pelican Coil Remover	50
400318	Pelican Coil Remover	5



Pelican Thread Retriever

The Pelican Thread Retriever is designed to quickly and easily locate lost IUD threads and can also be used to sound the uterus for the presence of the intra-uterine device.

Seven specifically designed notches at 10mm intervals (on either side), combined with a long handle, allow for greater control and optimum retrieval of IUD threads.

Features

- ✓ Thumb recess for enhanced grip and control
- ✓ Flexible stem with smooth rounded tip to sound the uterus
- ✓ Seven 10mm interval notches for optimum retrieval of IUD threads
- ✓ Long handle (245mm) for ease of use



Manufactured
in the UK

Code	Description	Units
400565	Pelican Thread Retriever	25
400560	Pelican Thread Retriever	5

Cervical Dilator & Uterine Sound





Pelican Uterine Sound-Dilator

NEW

The Uterine Sound-Dilator (“USD”) is a strong but flexible, single-use, double-ended combined device designed for gentle dilation of the cervix and to assess the depth of the uterus.

The smooth and rounded ends are specially designed to allow gradual dilation of the cervix, whilst reducing risk of perforation.

This ‘2-in-1’ instrument, provides added convenience when fitting intra-uterine devices

Features

- ✓ Strong but flexible moulded plastic
- ✓ Ideal for coil-fitting procedures
- ✓ Smooth, rounded tips for reduced risk of perforation
- ✓ Clear and precise 1cm depth indicators
- ✓ 27.7cm long with a maximum diameter of 5.5mm




2 in 1



Manufactured in the UK

Code	Description	Units
401000	Uterine Sound-Dilator	50
401010	Uterine Sound-Dilator	10



Pelican Dilator

The Pelican Dilator is designed to gently dilate the cervical 'os' to facilitate access to the uterus, and is ideal for preparing the cervix for the introduction of an intra-uterine device.

The Pelican Dilator is tapered to allow gradual dilation in 1mm increments from 2mm to 6mm (using two devices).

Features

- ✓ Strong but flexible moulded plastic
- ✓ Central thumb hold for ease of use
- ✓ Tapered ends allow gradual cervical dilation
- ✓ Smooth, rounded tips for reduced risk of perforation



Manufactured
in the UK

Code	Description	Units
400342	Pelican Dilator 3/4mm	100
400343	Pelican Dilator 3/4mm	10
400344	Pelican Dilator 5/6mm	100
400345	Pelican Dilator 5/6mm	10
400347	Pelican Dilator 7/8mm	100
400348	Pelican Dilator 7/8mm	10



Pelican Uterine Sound

The Pelican Uterine Sound is specifically designed to confirm uterine depth prior to the insertion of an intra-uterine device. The rounded end ensures a gentle introduction of the device through the cervical 'os' into the uterine cavity.

The Pelican Uterine Sound is a slender, flexible gauge with 1cm indicators to allow immediate visual confirmation of uterine depth.

Features

- ✓ Notches allow measurement of the uterine depth
- ✓ White in colour to aid visualisation
- ✓ Mirrors the metal equivalent for ease of use
- ✓ Increased patient comfort due to design
- ✓ Unique product identification number for traceability



Manufactured
in the UK

Code	Description	Units
400310	Pelican Uterine Sound	100
400311	Pelican Uterine Sound	10

Obstetrics Range





Pelican FETOfit®

NEW

FETOfit® is a disposable device specifically designed to secure both 'TOCO' and 'CTG' transducers quickly and easily for strapless foetal monitoring.

Unlike standard circumferential straps, FETOfit's hypoallergenic adhesive pads adhere directly to the front of the abdomen.

The flexible material moulds to body contours, offering a secure fit whilst promoting maternal comfort and mobility whilst facilitating access to the back in the event of an epidural.

Features

- ✓ Universal attachment to suit all transducers
- ✓ Quick and easy to re-position
- ✓ Easy access for epidural
- ✓ Disposable device for single-person use
- ✓ Secure fit
- ✓ Safeguards patient comfort & mobility



Manufactured in the UK

Code	Description	Units
400625	FETOfit®	25



PELIhook Amniotic Membrane Perforator “Amnihook”

The PELIhook Amniotic Membrane Perforator is designed to safely rupture the amniotic membrane whilst minimising the risk of injury to the foetus.

The sharp hook cleanly snags the membrane, and the rounded tip ensures that the hook does not inadvertently come into contact with maternal tissue.

Features

- ✓ The curved PELIhook is anatomically shaped to allow ease of use
- ✓ Designed to safely rupture the membrane
- ✓ Textured thumb grip on the handle
- ✓ Option of a straight PELIhook for optimal control and comfort
- ✓ Sterile and single-use
- ✓ Latex-free



Manufactured
in the UK

Code	Description	Units
SH2	PELIhook Amniotic Membrane Perforator - Straight	100
KH2	PELIhook Amniotic Membrane Perforator - Curved	100

Model-Med: Anatomical Training Models



Pelican Feminine Healthcare are the authorised distributors of the renowned Model-Med range of anatomical models for Gynaecology and Obstetrics medical simulations.


All models are made of superior materials that offer a far more realistic simulation of human tissue than other medical models available. The soft synthetic flesh will not rip or tear. The models do not require any special lubricants; any standard water based lubricant, such as PELIjelly, can be used.

- ✓ **Durable**
made of a very soft, highly elastic, and very durable material.
- ✓ **Easy to clean**
simple maintenance instructions ensures the models will last a long time.
- ✓ **Easy to repair**
easy to repair if accidentally damaged.

The range includes models for:

- Assessment of cervical dilatation and foetal presentation
- Full birth simulation
- Instrumental delivery of the foetus
- Full gynaecological model for speculum and bi-manual examination

The Model-med range is supplied with a one year warranty.



'Sophie'

Full Birth Simulator

Sophie's Mum is designed for the simulation of cephalic birth, breech birth, and shoulder dystocia. Sophie's Mum has a full pelvis, the joints of which are slightly loosened to simulate the movement in those joints during childbirth. The baby must travel around the 'Curve of Carras' and will be born in an accurate simulation of the birth.

Sophie is a new-born baby weighing 4.0 kg (8.8lbs) and measuring 45cm (18 inches) in length. She contains a fully articulated metal skeleton that behaves the same as a human skeleton.



Full care and handling instructions can be downloaded online from the Pelican Feminine Healthcare website.

Features

- ✓ Provides a highly realistic birthing experience
- ✓ Sophie includes a fully articulated metal skeleton
- ✓ Sophie's mum simulates the movement of that in childbirth



Code	Description	Units
OT002	Full Birth Base Unit	1
FT001	Foetus (Sophie)	1
UT008	Post Partum Uterus	1
PL001	Placenta	1
UT009	Post Partum Haemorrhage Uterus	1



'Charlie'

Obstetric Trainer

The Charlie Obstetric Trainer consists of Charlie (the head), Charlie's mum the base unit), and three different cervix: primipara, multipara, and an oedematous cervix.

Charlie and his mum are made from a remarkably soft synthetic flesh that provides a highly realistic experience for obstetric training.



Full care and handling instructions can be downloaded online from the Pelican Feminine Healthcare website.

Features

- ✓ Provides a highly realistic experience
- ✓ Available with 3 different types of cervix
- ✓ Both models are made from a soft synthetic flesh material
- ✓ Tough flesh material does not tear or split

Code	Description	Units
OT001	Charlie's Mum Obstetric Trainer Base Unit	1
FH001	Foetal Head (Charlie)	1
CM001	Cervix Multipara	1
CP001	Cervix Primipara	1
CO001	Cervix Oedematous	1



Lucy & Lucy's Mum Instrumental Delivery Trainer

Lucy's Mum is specially designed to practice instrumental delivery, particularly vacuum assisted delivery.

Lucy's head is the average shape of the foetal head. It can be used to practice forceps delivery, but is particularly suited to vacuum delivery. Together with the Lucy's Mum model it provides an excellent training platform for medical professionals.




Full care and handling instructions can be downloaded online from the Pelican Feminine Healthcare website

Features

- ✓ Provides a highly realistic experience
- ✓ Lucy's head mimics the shape and weight of an actual foetal head
- ✓ Can be used with forceps
- ✓ Ideal for practising vacuum assisted delivery
- ✓ Other skin tone options available



Code	Description	Units
OT004	Lucy's Mum Instrumental Delivery Base Unit	1
FT002	Lucy Vacuum Extraction Head	1
LMUK001	Lucy's Mum Instrumental Delivery Base Unit & Head (Light skin tone)	1
LMUK002	Lucy's Mum Instrumental Delivery Base Unit & Head (Dark skin tone)	1



Gynaecological Trainer

The Model-Med Gynaecological Trainer allows a totally convincing manual and visual examination of the female reproductive system. A speculum can be used on the model to simulate procedures, such as a pap smear.

Using a unique clip-fit system, different uteri can be swapped in and out of the base unit. Seven types of uterus are available.



Full care and handling instructions can be downloaded online from the Pelican Feminine Healthcare website.

Features

- ✓ Provides a highly realistic experience
- ✓ Available with 7 different types of uterus
- ✓ Suitable for use with a speculum for practising smear taking
- ✓ Synthetic flesh material is soft but tough and will not tear or split



Code	Description	Units
GN001	Trainer Base Model	1
UT001	Small Uterus	1
UT002	Medium Uterus	1
UT003	Large Uterus	1
UT004	Fibroid Uterus	1
UT005	Retroverted Uterus	1
UT006	Hollow Uterus	1
UT007	Ovarian Cyst Uterus	1
CX001	Extra Clip	1

NEW

FETOFit®

Strapless foetal monitoring

FETOFit® is a disposable device specifically designed to secure both 'TOCO' and 'CTG' transducers quickly and easily for foetal monitoring.

Unlike standard circumferential straps, FETOFit's **hypoallergenic** adhesive pads adhere directly to the front of the abdomen for a **secure fit** with patient **comfort** and **mobility**.



Order information

Description:	FetoFit®
Units:	25
Product Code:	400625
NHS Code:	FGR450

Pelican Feminine Healthcare Ltd.
 Greypoint, Cardiff Business Park, Cardiff CF14 5WF
0800 013 1220 or contactus@pelicanfh.co.uk
www.pelicanfh.co.uk

- ✓ Secure fit
- ✓ Quick & easy to re-position
- ✓ Easy access for epidural
- ✓ Safeguards patient comfort & mobility


www.pelicanfh.co.uk

Pelican Feminine Healthcare are delighted to announce the launch of their brand new website. Browsing across mobile, desktop and tablet devices has never been so easy!

A commitment to delivering market-leading Gynaecology and Obstetrics products with excellent service to all our customers, at **www.pelicanfh.co.uk** you can find:

- ✓ New clean and simple design makes navigating even easier
- ✓ Quickly access the Pelican Feminine Healthcare's full product range
- ✓ High resolution product image gallery
- ✓ Latest product news, events and company blog
- ✓ Even more downloadable content
- ✓ Easy sign-up for latest offers and product updates

Follow us on Twitter and LinkedIn today for all the latest news and offers!



Proud to supply market
leading Gynaecology
and Obstetrics products
to health professionals
worldwide

Freephone: 0800 013 1220

E-mail: contactus@pelicanfh.co.uk

www.pelicanfh.co.uk

Pelican Feminine Healthcare





Pelican Feminine Healthcare Ltd.

Greypoint, Cardiff Business Park, Cardiff CF14 5WF

0800 013 1220 contactus@pelicanfh.co.uk

www.pelicanfh.co.uk