



# High Volume

## HEMORRHOIDAL CIRCULAR STAPLER

Specifically studied for the treatment  
of hemorrhoidal prolapse and  
for the S. T. A. R. R. technique





Patent design:  
Transparent cartridge,  
enhanced visibility and  
control over captured tissue



2 big casing ports  
offer extra visibility  
when pulling tissue  
into the cartridge

Anvil



The biggest case  
in the market for  
a 34 hemorrhoidal  
stapler: 27cc



Markings, indicates  
the distance  
of anastomotic stoma  
from the anal verge



Obturator



Purse-string  
Suture Anoscope



Suture  
Threader



Circular Anal  
Dilator

# HIGH VOLUME HEMORRHOIDAL CIRCULAR STAPLER



Enlarged indicator window, convenient for precise adjustment of staples closure height



Closing knob: easy to open and close, convenient for operation



External safety

Sound feedback, indicates the stapler is ready to be fired, safer operation

Firing handle

Patent design: Dual safety, avoids misoperation, much safer



Butterfly Anal Dilator

# HIGH VOLUME HEMORRHOIDAL CIRCULAR STAPLER



Product	Code	Description
	FCSSWBE31	Disposable Hemorrhoidal Circular Stapler 31mm, dual safety
	FCSSWBE32	Disposable Hemorrhoidal Circular Stapler 32mm, dual safety
	FCSSWBE33	Disposable Hemorrhoidal Circular Stapler 33mm, dual safety
	FCSSWBE34	Disposable Hemorrhoidal Circular Stapler 34mm, dual safety
	FCSSWAE31	Disposable Hemorrhoidal Circular Stapler 31mm, dual safety, with transparent cartridge
	FCSSWAE32	Disposable Hemorrhoidal Circular Stapler 32mm, dual safety, with transparent cartridge
	FCSSWAE33	Disposable Hemorrhoidal Circular Stapler 33mm, dual safety, with transparent cartridge
	FCSSWAE34	Disposable Hemorrhoidal Circular Stapler 34mm, dual safety, with transparent cartridge
	FCSSME32	High Volume Disposable Hemorrhoidal Circular Stapler 32mm, dual safety, With transparent cartridge and sound feedback
	FCSSME34	High Volume Disposable Hemorrhoidal Circular Stapler 34mm, dual safety, with transparent cartridge and sound feedback